

[برگشت \(https://eqclد.com/IACLD/Panel/Labs/Answers/34409\)](https://eqclد.com/IACLD/Panel/Labs/Answers/34409)

چاپ

دوره چهل و هفتم زمستان 1402

### B3 Results

#### نمونه مجهول پروتئین ادرار

کد نمونه 2487	EQAP-10814	آزمایشگاه نوبل
---------------	------------	----------------

عملکرد عالی (رنگ آبی)	DI	< 0.5
عملکرد مطلوب (رنگ سبز)	DI	0.5 - 1.0
نتیجه قابل قبول اما در حدود مرز (رنگ زرد)	DI	1.0 - 2.0
نیاز به بررسی روش و کنترل کالیبراسیون (رنگ نارنجی)	DI	2.0 - 3.0
وجودنقص و نیاز به پیگیری فوری (رنگ قرمز)	DI	> 3.0
نتیجه گزارش شده خارج از محدوده است (Out Of Range)	DI	> 10.0

برای مشاهده آنالیز آماری تمامی همگروه ها، در جدول زیر بر روی نام آنالیت مورد نظر کلیک کنید.

Analyte	Unit	Lab Result	PeerGroup Count	PeerGroup Mean	PeerGroup %CV	Lab DI
<b>Urine Protein</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183771">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183771</a> )	mg/L	364.00	237	368.73	9.67	<b>0.13</b>
<b>Protein/Cr Ratio</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183802">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183802</a> )	mg/g Cr	539.00				
<b>Protein (24h)</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183824">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183824</a> )	mg/24h	728.00				
<b>Creatinine B3</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183834">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183834</a> )	mg/dL	67.50	7	66.00	23.68	<b>0.23</b>
<b>Creatinine 24</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183868">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183868</a> )	mg/24h	0.07				
<b>Calcium</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183877">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183877</a> )	mg/dL	5.00	8	5.36	40.23	<b>0.68</b>
<b>Calcium 24</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183908">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183908</a> )	mg/24h	100.00				
<b>Calcium/Creatinine Ratio</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183927">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183927</a> )	mg/mg	0.07				
<b>Uric Acid B3</b> ( <a href="https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183945">https://eqclد.com/IACLD/Panel/Analysis/Report/34409/595/183945</a> )	mg/dL	14.20				

<b>Uric Acid 24</b> ( <a href="https://eqclid.com/IACLD/Panel/Analysis/Report/34409/595/183962">https://eqclid.com/IACLD/Panel/Analysis/Report/34409/595/183962</a> )	mg/ 24h	284.00				
<b>UricAcid/Creatinine Ratio</b> ( <a href="https://eqclid.com/IACLD/Panel/Analysis/Report/34409/595/183978">https://eqclid.com/IACLD/Panel/Analysis/Report/34409/595/183978</a> )	mg/ mg	0.21				

---

Copyright © 2024 IACLD. All rights reserved.